

**1.5T ADNI GE 12.0 M4 Software, QC Phantom  
TwinSpeed Gradient and 8-channel Brain Coil**

---

Set up the phantom according to the ADNI instruction manual. Use the 8-channel brain coil.

Use the electronic ADNI protocol that has been loaded on your scanner.

**NOTE: The "adni" pulse sequence used for series 2 is a works in progress pulse sequence that ADNI sites obtain from GE Healthcare.**

<b>SERIES</b>	1. 3 plane loc.	scan plane mode SAT	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	coil 8hrbrain		3-plane	matrix/nex	256 / 128 / 1	
	etl		(Whole Body gradient)	fov (cm)	26	
	scan time :13			slice/space	5 / 5	
			autoshim	On		
comments	Use 8-channel brain coil.					

<b>SERIES</b>	2. Sag MPRAGE	scan plane mode pulse seq image opts. psd name	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>		
	coil 8hrbrain		Sag	matrix/nex	192 / 192 / 1		
	#echos te		SCAN TIMING	3D ( ZOOM gradient)	phase fov		
				SPGR	locs/pause		
Prep time	1	EDR, IrPrep, Fast	freq. direct.	S/I			
flip angle	min full	adni	fc direct				
etl	1000	User CVs	<i>ADDITIONAL PARAMETERS</i>		phase corr		
bw1/bw2	8		Image acq. delay = 0	autoshim	On		
scan time	15.63		turbo mode = 1	fov	<i>SCANNING RANGE</i>		
	7:42		swap fat chem shift = 1	slice/space			24
		slice resolution = 100%		1.3mm 170 locs/slab			
		mprage mode = 1					
comments		TD or mprageTR=2400					
	Cover the phantom						